

## Non-invasive Mask for Hospital

**13521:** 834 NON VENTED+AAV FULL FACE MASK, S

**13522:** 834 NON VENTED+AAV FULL FACE MASK, M

**13523:** 834 NON VENTED+AAV FULL FACE MASK, L

\*Carton pack: 1pc / pack



\* Flexible, one-piece forehead pad evenly conforms to the forehead contour



\* Innovative adjustable forehead support provides stability with easy one-handed adjustment



\* Medical grade liquid silicone full face cushion fits just right and firmly on the face

- Quick assembling design offers user more convenience when taking off / assembling the cushion
- Sampling port allows measuring pressure or ETCO2
- Minimum components will lower the risk of losing components and easier for cleaning
- Push release / attach buckle, easier for user to take off mask
- 360 degree rotational ball in elbow



Vented+AAV



Non-vented



Non-vented+AAV

\* Optional Elbow Connectors

Manufactured by  
HSINER CO., Ltd.

No.312, Jhongshan Rd., Shengang Dist., Taichung City 429, Taiwan

Tel: 886-4-25152480

Fax: 886-4-25152482

Web Site: <http://www.hsiner.com>

