


Néa, the only foam-free CPAP

Néa optimises patient safety and comfort thanks to its patented Sefam features




And for the record, Néa is the size of this document.



Head office: Sefam, 144 avenue Charles de Gaulle - 92200 Neuilly-su-Seine - FRANCE
Manufacture: Sefam, 10 allée Pelletier Doisy - 54600 Villers-Lès-Nancy - FRANCE





BR0231-01 / 06-2024



A simple access to essential functions



Full-size screens



Circuit Select
Estimation and consideration of the resistance of the patient's circuit



Mask Select
Taking into account the intentional leak of the used mask



Mask Off
Recognition of disconnected mask



Connectivity
BLE, GSM, WiFi



iRamp
Intelligent ramp



Intelligent Start
Automatic start



CC+
Comfort mode



more information



Néo is a Continuous Positive Airway Pressure (CPAP Class IIa) medical device connected by a flexible tube to a mask worn by patients suffering from obstructive sleep apneas and hypopneas syndrome (OSAHS) and weighing over 30kg. The device delivers a continuous flow of air to avoid upper airway obstruction during sleep. It is a regulated healthcare product and, as such, bears the CE mark. Please refer to the device's user manual for information on recommended use, contraindications, side effects and correct, safe and effective use.

Néo is compatible with Sefam Connect, Sefam Analyze, Sefam Access Lite and Pro