

Geriatric Chair

Your Best Choice

Full Recline and with 3 Position Locking System

Item Code: UM-R7001



Features:

- full recline from back angle 104° - 165° and with safety locking in 104°, 143° and 157°
- with 3 position safety locking system for recline: upright, deep recline and elevated legrest/ footrest
- safety locking system secures the chair in the desired position
- 125mm castors make it easy to move and secure the recliner
- CMHR high density foam for the seat and back pad
- neutral colored, sturdy steel frame with non-toxic epoxy coating, and easy to clean and maintain
- blow mould meal tray locks in 5 positions
- two side with urine bag hanger
- PVC coating pedal with sliding footrest

Optional Accessories:

- 50mm heavy duty nylon safety belt
- 165mm height pommel
- 4 hooks height adjustable drip stand
- 'U' shape headrest
- contour lumbar cushion with velcro strap



Upright Position



Elevated Legrest



Deep Recline



Free Recline Position



3 Position Locking system



'U' shape Urine Bag Hanger

Specifications:

- Overall Depth: 920 - 1700mm
- Overall Height: 1170 - 770mm
- Overall Width: 650mm
- Seat Width: 480mm
- Seat Depth: 450mm
- Seat Pad Width: 460mm
- Seat height (Upright Position): 460mm in front
- Back Height (measure from seat level): 900mm
- Back Pad Height: 850mm
- Back Pad Width: 460mm
- Armrest Pad: 760mm (L) x 90mm(W)
- Flame standard: BS7176:2007 medium hazard
- Upholstery: Vinyl
- Loading: 120 kgs